



MERCK

THE SAFETY SENTINELS

Environmental Monitoring Solutions
for Pharmaceutical Industry



The life science business of Merck operates
as MilliporeSigma in the U.S. and Canada.

Millipore®

Preparation, Separation,
Filtration & Monitoring Products

THE SAFETY SENTINELS

MERCK



INTRO



APPLICATIONS



CLEANROOM
GRADES



SERVICE



ONE

STEP AHEAD!

in Environmental Monitoring Solutions

CONTACT



THE SAFETY SENTINELS



INTRO



APPLICATIONS



CLEANROOM
GRADES



SERVICE

Surface & Personnel Monitoring

Passive Air Monitoring

Active Viable Air Monitoring

Compressed Gas Monitoring

Media Fill / Aseptic Process Simulation



ONE

STEP AHEAD!

in Environmental Monitoring Solutions

CONTACT



THE SAFETY SENTINELS

MERCK



INTRO



APPLICATIONS



CLEANROOM
GRADES



SERVICE

Cleanroom Grade A/B

Cleanroom Grade C/D

Isolator / RABS



ONE
STEP AHEAD!

in Environmental Monitoring Solutions

CONTACT



ON GUARD



GMP compliant solutions for process monitoring in pharmaceutical manufacturing

Safety matters. The manufacture of pharmaceutical products is performed under controlled conditions. Microbial monitoring is an important part of proving the manufacturing process is under control, especially in aseptic production. The risk-based implementation of sampling plans as well as subsequent trend analysis of monitoring results, this helps to detect deviations from qualified status, prepare root cause analysis and follow up with appropriate CAPAs.



Learn more on
our Website



100%

OUR SERVICES FOR YOU

Benefit from our all-around services offer from validation and training to maintenance and calibration.

- Validation package: validation protocols and on-site IQ/OQ service
- Service plans: preventative & curative maintenance
- Traceability & calibration accuracy

Learn more on
our Website





ONE


STEP AHEAD!

in Environmental Monitoring Solutions

Learn more on
our Website



ACCURACY



setting standards



MINIMIZES RISK



Accuracy



Consistent and reliable high performance

Our products are manufactured for high quality assured raw materials, each raw material is from a proven and audited supplier. The use of standardized manufacturing practices ensures high quality products. These conditions ensure a high quality and reliable product that makes your result trending data accurate, reliable and consistent.



setting standards

for ease of use and safety

Our expertise and constant quality improvement programs coupled with a high level of validation and scientific data advances the environmental monitoring products to a world leading position.



Minimizes Risk



to make your work more secure

Our product range is designed to minimize the risk of cross-contamination as well as disruption of unidirectional air flow during active air monitoring up to the highest level of cleanroom grades and isolators.



SURFACE & PERSONNEL MONITORING



**HYCON®
Contact Slides**

**Cleanroom Grade
A/B**

**Cleanroom Grade
C/D**



**ICR/ICR+
Contact Plates**
55 mm

**Cleanroom Grade
A/B**

Isolator / RABS



ICR Swabs

**Cleanroom Grade
A/B**

Isolator / RABS



**IsoBag®
Rapid Transfer Bag**

**Passive Air
Monitoring**

Isolator / RABS



**RT/RT+
Contact Plates**
55 mm

**Cleanroom Grade
C/D**

HYCON® Contact Slides

for monitoring curved surfaces and personnel



Description

Benefits

Ordering Information

HYCON® Contact Slides are designed for monitoring flat and curved surfaces as well as personnel. They are available in standard single packaging or in a double-bagged, gamma-irradiated format.

Both the standard and the gamma-irradiated **HYCON® Contact Slides** are rectangular in shape and provide a surface of 25 cm². Individual sealing of the primary packaging allows opening single contact slides at a time to reduce waste.

The slide's body is slightly flexible to allow efficient sampling of curved surfaces. All slides are equipped with a cover slide for closure after use, increasing safety during both transport and incubation.



[Link to our Website](#)



[Download Brochure](#)



HYCON® Contact Slides

for monitoring curved surfaces and personnel



Description

Benefits

Ordering Information

- For critical or non-critical environments
- Designed for use on flat and curved surfaces
- Safe transport and incubation



HYCON® Contact Slides

for monitoring curved surfaces and personnel



Description

Benefits

Ordering Information

Product	Description	Content	Ord. No.
Contact Slides Blank Strip	Kit for use with HYCON®	20	1.44024.0020
HYCON® Contact Slides	Gamma-irradiated, double-bagged		
Contact Slides DE-y	Gamma-irradiated, modified Dey Engley Agar (D/E Agar), double-bagged; for determination of total count in aseptic environments; with neutralizers (lecithin, sorbitol monooleate complex, sodium thioglycolate, sodiumbisulfite) and growth supplements.	20	1.44021.0020
Contact Slides SDX-y	Gamma-irradiated Sabouraud Dextrose Agar with modified pharmacopoeia formulation, double-bagged; for determination of yeasts and molds in aseptic environments; with neutralizers (Tween® 80, L-a-phosphatidylcholine) and growth supplements.	20	1.44016.0020
Contact Slides TC-y	Gamma-irradiated Tryptic Soy Agar, double-bagged; for determination of total count in aseptic environment; with neutralizers (lecithine, sorbitol monooleate complex, Tween® 80) and growth supplements.	20 100	1.44022.0020 1.44022.0100



HYCON® Contact Slides

for monitoring curved surfaces and personnel



Description

Benefits

Ordering Information

Product	Description	Content	Ord. No.
HYCON® Contact Slides	Standard, single-bagged		
Contact Slides C	Mac Conkey Agar; for determination of coliform bacteria; with growth supplements.	20	1.44019.0020
Contact Slides SDX	Sabouraud Dextrose Agar with modified pharmacopoeia formulation; for determination of yeasts and molds; with neutralizers (Tween® 80, L-a-phosphatidylcholine) and growth supplements.	20	1.44017.0020
Contact Slides TC	Tryptic Soy Agar; for determination of the total count; with neutralizers (lecithin, L-histidine, sorbitol monooleate complex, Tween® 80) and growth supplements.	20 100	1.44023.0020 1.44023.0100
Contact Slides YM	Rose Bengal Agar with streptomycin and chloramphenicol; for determination of yeasts and molds; with growth supplements.	20	1.44018.0020



ICR/ICR+ Contact Plates

for isolators and critical cleanrooms



Description

Benefits

Related Products

Further Reading

Accessories

Ordering Information

The **ICR/ICRplus series of contact plates** provide maximum safety for use in isolators and cleanrooms.

- Aseptically filled and irradiated in final packaging to minimize contamination.
- Full traceability due to Data Matrix barcode on each plate.
- Available with standard TSA/SDA formulations with wide variety of neutralizers for either disinfectants or β -lactam antibiotics.

ICRplus products:

- Smart "CLOSED-VENT" locking system allowing safe transport and optimized incubation under anaerobic conditions.



Link to our Website



Download Brochure



Request a demo



ICR/ICR+ Contact Plates

for isolators and critical cleanrooms



Description

Benefits


Related Products

Further Reading

Accessories

Ordering Information

**CLOSED- or
VENT-Closure**



Maximum safety

- Filled under aseptic conditions, gamma-irradiated in final packaging
- Safe transport thanks to H₂O₂-impermeable triple-sleeve packaging
- Data Matrix barcode to ensure full traceability of each individual plate

Maximum convenience

- Storage at room temperature to store at site of use
- Long shelf life to reduce incoming goods controls
- SDA in colored plates to minimize risk of mix-ups

Additional benefits for ICRplus Contact Plates

- Safe transport through integrated locking systems
- CLOSED-VENT system for different incubation conditions

ICR/ICR+ Contact Plates

for isolators and critical cleanrooms



Description

Benefits

Related Products

Further Reading

Accessories

Ordering Information

ICR/ICRplus Settle Plates

Our gamma-irradiated and triple-bagged settle plates for monitoring of isolators and critical cleanrooms.



Link to our Website



ICR/ICR+ Contact Plates

for isolators and critical cleanrooms



Description

Benefits

Related Products

Further Reading

Accessories

Ordering Information

ICR Swabs

Easy presence-absence test on difficult to access surfaces.



Link to our Website



ICR/ICR+ Contact Plates

for isolators and critical cleanrooms



Description

Benefits

Related
Products

Further
Reading

Accessories

Ordering
Information



Additional information

- AN4476 – Neutralization Efficacy of ICR Plates
- Reworked Anios AN
- SDA AN
- AN2005 – Anaerobic Monitoring using our lockable TSA+LTHTh ICR+ plates

ICR/ICR+ Contact Plates

for isolators and critical cleanrooms



Description

Benefits

Related Products

Further Reading

Accessories

Ordering Information

Product

**Plate Rack
for Contact Plates**

Description

Stainless Steel Rack
for 10 x 55 mm plates

Ord. No.

1.46502.0001



ICR/ICR+ Contact Plates

for isolators and critical cleanrooms



Description

Benefits

Related Products

Further Reading

Accessories

Ordering Information

Product	Description	Content	Ord. No.
ICR Contact Plates	55 mm , triple-bagged, gamma-irradiated, non-lockable		
SDA Contact +LTHTh – ICR	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate; filled in pink plates	20	1.46201.0020
		200	1.46201.0200
TSA Contact +LT – ICR	with neutralizers lecithin, Tween® 80	20	1.46195.0020
		200	1.46195.0200
TSA Contact +LTHTh – ICR	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate	20	1.46231.0020
		200	1.46231.0200
TSA LTHThio cont. – ICR	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate	20	1.46797.0020
		200	1.46797.0200



ICR/ICR+ Contact Plates

for isolators and critical cleanrooms



Description

Benefits

Related Products

Further Reading

Accessories

Ordering Information

Product	Description	Content	Ord. No.
ICRplus Contact Plates	55 mm , triple-bagged, gamma-irradiated, lockable		
Chocolate Contact Agar +LTH – ICR+	with neutralizers lecithin, Tween® 80, histidine, supports growth of fastidious aerobic and anaerobic microorganisms	200	1.46555.0020
SDA Contact +LTHTh – ICR+	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate; filled in pink plates	20 200	1.46501.0020 1.46501.0200
SDA Contact +LTHTh selective – ICR+	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate, irradiation-resistant antibiotics to inhibit bacterial growth; filled in pink plates	20 200	1.46538.0020 1.46538.0200
TSA Contact +LT – ICR+	with neutralizers lecithin, Tween® 80	20 200	1.46552.0020 1.46552.0200
TSA Contact +LT+ Cephase – ICR+	with neutralizers lecithin, Tween® 80, β-lacatamses for inactivation of penicillins, all generations of cephalosporins and carbapenems	200	1.46539.0200

ICR/ICR+ Contact Plates

for isolators and critical cleanrooms



Description

Benefits

Related Products

Further Reading

Accessories

Ordering Information

Product	Description	Content	Ord. No.
ICRplus Contact Plates	55 mm , triple-bagged, gamma-irradiated, lockable		
TSA with LTHThio cont. – ICR+	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate	20	1.46783.0020
		200	1.46783.0200
TSA with Penase LTHThio cont. – ICR+	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate, penase for inactivation of penicillins	20	1.46798.0020
		200	1.46798.0200
Vegetable pepton agar LTHThio cont. – ICR+	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate	20	1.46803.0020
		200	1.46803.0200
Neutralizer A – Contact – ICR+	neutralizers for increased spectrum of sanitizers	20	1.46697.0020
		200	1.46697.0200

ICR Swabs



for surface and personnel monitoring in isolators and cleanrooms

Description

Benefits

Related Products

Accessories

Ordering Information

Our **ICR swabs** are designed for absence-presence testing on dry, difficult-to-access surfaces and for personnel monitoring. The patented snap-valve design of the swab makes handling both convenient and safe. The low-abrasion swab tip material minimizes the release of particles during use. The tip's moistening solution does not leave any growth promoting residues on the sampled surfaces.

The swab's TSB medium is supplemented with neutralizers to ensure optimal growth conditions, even in the presence of VHP residues in isolators.

A stainless-steel rack helps to safely and easily transport larger numbers of ICR Swabs. It holds 36 tubes and is sized 196 x 99 x 104 mm.



[Link to our Website](#)



[Download Brochure](#)



[Request a demo](#)



ICR Swabs



for surface and personnel monitoring in isolators and cleanrooms

Description

Benefits

Related Products

Accessories

Ordering Information



- Suitable for cleanroom and isolator use
- Reaches difficult to access areas
- All-in one design for low contamination risk and easy handling
- Residue-free sampling
- Qualitative microbial presence/absence test

ICR Swabs

for surface and personnel monitoring in isolators and cleanrooms



Description

Benefits

Related Products

Accessories

Ordering Information

ICR/ICRplus Contact Plates

Gamma-irradiated and triple-bagged contact plates are suitable for use in critical cleanrooms, RABS and isolators.



Link to our Website



ICR Swabs

for surface and personnel monitoring in isolators and cleanrooms



Description

Benefits

Related Products

Accessories

Ordering Information

Product


Description

Ord. No.

Swab Rack

Stainless steel

1.46530.0001



ICR Swabs



for surface and personnel monitoring in isolators and cleanrooms

Description

Benefits

Related Products

Accessories

Ordering Information

Product	Content	Ord. No.
ICR Swabs	100	1.46529.0100



IsoBag®

for rapid transfer of culture media plates into isolator



Description

Benefits

Properties

Related Products

Ordering Information

Production workflows are always at their most efficient when all tools are ready and available. The **IsoBag**® rapid transportation and packaging bag offers a reliable and easy way to have environmental monitoring plates right at hand in the isolator.

IsoBag® speeds up environmental monitoring in aseptic production isolators as it provides ready-to-use gamma-irradiated contact or settle plates for immediate use; the bag is simply mounted to the 190 mm alpha port of the isolator, and the required amount of plates can be transferred into the isolator for immediate use.

Saving time spent on decontamination increases operational productivity. In addition, no valuable space is wasted in the isolator to store the required amount of plates between decontamination cycles.



[Link to our Website](#)



[Download Brochure](#)



[Request a demo](#)



IsoBag®

for rapid transfer of culture media plates into isolator



Description


Benefits

Properties

Related Products

Ordering Information

- **No need to decontaminate plates' primary bags**
- **Containment** is maintained throughout use
- Plates are held **ready for immediate use**
- **No storage space** needed within isolator



IsoBag®

for rapid transfer of culture media plates into isolator



Description

Benefits

Properties

Related
Products

Ordering
Information



- Pre-packed DPTE BetaBag® with ready-to-use single-bagged contact or settle plates
- Beta port connectable to any 190 mm DPTE® Alpha port
- Finished product gamma-irradiated for immediate use
- Contained plates traceable to corresponding batch of ICR plates
- Available with our TSA + LTHTh in lockable and non-lockable contact and settle plates
- Multiple connections possible
- Sterile bags for safe transport of used plates included

IsoBag®

for rapid transfer of culture media plates into isolator



Description

Benefits

Properties

Related Products

Ordering Information

ICR/ICRplus Contact Plates

Gamma-irradiated and triple-bagged contact plates are suitable for use in critical cleanrooms, RABS and isolators.



Link to our Website



IsoBag®

for rapid transfer of culture media plates into isolator



Description

Benefits

Properties

Related Products

Ordering Information

ICR/ICRplus Settle Plates

Our gamma-irradiated and triple-bagged settle plates for monitoring of isolators and critical cleanrooms.



Link to our Website



IsoBag®

for rapid transfer of culture media plates into isolator



Description

Benefits

Properties

Related Products

Ordering Information



ICR Swabs

Easy presence-absence test on difficult to access surfaces.



Link to our Website



IsoBag®

for rapid transfer of culture media plates into isolator



Description

Benefits

Properties

Related Products

Ordering Information

Product	Description	Ord. No.
IsoBag® with Contact Plates	for 190 mm alpha-port	
IsoBag® TSA+LT Contact	10 x 10 Contact Plates 146195, 55 mm irradiated	1.46784.0100
IsoBag® TSA+LTHTh Contact	55 mm Contact Plates 146231, irradiated	1.46754.0100
IsoBag® TSA+LTHThio Contact+	55 mm lockable Contact Plates 146783, irradiated	1.46753.0100
IsoBag® TSA+LT Contact+	55 mm lockable Contact Plates 146552, irradiated	1.46817.0100
IsoBag® with Settle Plates	for 190 mm alpha-port	
IsoBag® TSA+LT Settle+	8 x 10 lockable Settle Plates 146684, 90 mm, irradiated	1.46785.0080
IsoBag® TSA+LTHTh Settle	90 mm Settle Plates 146069, irradiated	1.46756.0080
IsoBag® TSA+LTHTh Settle+	90 mm lockable Settle Plates 146683, irradiated	1.46755.0080
IsoBag® TSA Settle	90 mm Settle Plates 146001, irradiated	1.46814.0080

RT/RT+ Contact Plates



for microbial monitoring of air, surfaces and personnel

Description

Benefits

Ordering Information



Microbial monitoring of air, surfaces and personnel requires culture media that ensure reproducible detection results. RT tryptic soy agar (TSA) and sabouraud dextrose agar (SDA) reliably detect bacteria, yeasts and molds in non-critical areas.

The TSA and SDA **RT and RTplus contact plates** are suitable to monitor surfaces for bacterial and fungal contamination.

Selective **RT contact plates** support identification of Staphylococcus ssp., Escherichia coli, bile-tolerant Gram-negative bacteria, yeasts and molds.

Link to our Website



Download Brochure



Request a demo



RT/RT+ Contact Plates





for microbial monitoring of air, surfaces and personnel

Description

Benefits

Ordering Information

**CLOSED- or
VENT-Closure**



- **Controlled high-quality raw materials**
- **Data Matrix barcode** to ensure full traceability of each individual plate
- **Filled under aseptic conditions**
- **Long shelf life** to reduce incoming goods controls
- **Storage at room temperature** to keep at site of use

Additional benefits for RTplus Contact Plates

- Safe transportation thanks to **integrated locking systems**
- **CLOSED-VENT system** for different incubation conditions

RT/RT+ Contact Plates



for microbial monitoring of air, surfaces and personnel

Description

Benefits

Ordering Information

Product	Description	Content	Ord. No.
RT Contact Plates	55 mm , single-bagged, non-lockable		
Malt Extract Agar Contact – RT		20	1.46191.0020
Plate Count Agar Contact – RT		20	1.46154.0020
TSA Contact – RT		20	1.46240.0020
TSA Contact +LTH – RT	with neutralizers lecithin, Tween® 80, histidine	20	1.46200.0020
RTplus Contact Plates	55 mm , single-bagged, lockable		
SDA Contact+ Chloramphenicol – RT+	lockable plate	20	1.46549.0020
TSA Contact +LTH – RT+	with neutralizers lecithin, Tween® 80, histidine, lockable plate	20	1.46554.0020



PASSIVE AIR MONITORING



**ICR/ICR+
Settle Plates**

90 mm / 150 mm

**Cleanroom Grade
A/B**

Isolator / RABS



**IsoBag®
Rapid Transfer Bag**

**Surface & Personnel
Monitoring**

Isolator / RABS



**LI
Settle Plates**

90 mm

**Cleanroom Grade
C/D**

ICR/ICR+ Settle Plates

for isolators and critical cleanrooms



Description

Benefits

Related Products

Further Reading

Ordering Information

The **ICR/ICRplus series of settle plates** provide maximum safety for use in isolators and cleanrooms.

- Aseptically filled and irradiated in final packaging to minimize contamination.
- Full traceability due to Data Matrix barcode on each plate.
- 30 mL high filling volume to prevent desiccation under extended use.
- Available with standard TSA/SDA formulations with wide variety of neutralizers for either disinfectants or β -lactam antibiotics.

ICRplus products:

- Smart “CLOSED-VENT” locking system allowing safe transport and optimized incubation under anaerobic conditions.



Link to our Website



Download Brochure



Request a demo



ICR/ICR+ Settle Plates

for isolators and critical cleanrooms



Description


Benefits

Related Products

Further Reading

Ordering Information

CLOSED- or VENT-Closure



Maximum safety

- Filled under aseptic conditions, gamma-irradiated in final packaging
- Safe transport thanks to H₂O₂-impermeable triple-sleeve packaging
- Data Matrix barcode to ensure full traceability of each individual plate

Maximum convenience

- Storage at room temperature to store at site of use
- Long shelf life to reduce incoming goods controls
- SDA in colored plates to minimize risk of mix-ups
- High filling volume to allow prolonged exposure and incubation

Additional benefits for ICRplus Contact Plates

- Safe transport through integrated locking systems
- CLOSED-VENT system for different incubation conditions

ICR/ICR+ Settle Plates

for isolators and critical cleanrooms



Description

Benefits

Related Products

Further Reading

Ordering Information

ICR/ICRplus Contact Plates

Gamma-irradiated and triple-bagged contact plates are suitable for use in critical cleanrooms, RABS and isolators.



Link to our Website



ICR/ICR+ Settle Plates

for isolators and critical cleanrooms



Description

Benefits

Related Products

Further Reading

Ordering Information

MAS air sampler family

The industry standard for viable air sampling.



MAS-100 NT®



MAS-100 VF®



MAS-100 Iso MH®



MAS-100 Iso NT®



Link to our Website



ICR/ICR+ Settle Plates

for isolators and critical cleanrooms



Description

Benefits

Related
Products

Further
Reading

Ordering
Information



Additional information

- TN4909 – Air Sampling with pre-dried and fresh ICR agar plates
- AN2005 – Anaerobic Monitoring using our lockable TSA+LTHTh ICR+ plates
- SDA AN

ICR/ICR+ Settle Plates

for isolators and critical cleanrooms



Description

Benefits

Related Products

Further Reading

Ordering Information

Product	Description	Content	Ord. No.
ICR Settle Plates	90 mm , triple-bagged, gamma-irradiated, non-lockable		
SDA +LT – ICR 30 mL	with neutralizers lecithin, Tween® 80, filled in pink plates	120	1.46081.0120
SDA +LTHTh – ICR 30 mL	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate, filled in pink plates	20	1.46005.0020
		120	1.46005.0120
SDA selective +LTHTh – 30 mL ICR	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate, antibiotics against bacteria, filled in pink plates	120	1.46016.0120
TSA – ICR 30 mL		20	1.46001.0020
		120	1.46001.0120
TSA +LT – ICR 30 mL	with neutralizers lecithin, Tween® 80	20	1.46050.0020
		120	1.46050.0120
TSA +LTHTh – ICR 30 mL	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate	20	1.46069.0020
		120	1.46069.0120



ICR/ICR+ Settle Plates

for isolators and critical cleanrooms



Description

Benefits

Related Products

Further Reading

Ordering Information

Product	Description	Content	Ord. No.
ICR Settle Plates	90 mm , triple-bagged, gamma-irradiated, non-lockable		
TSA +LTHThio sedi. – ICR	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate	20 120	1.46786.0020 1.46786.0120
TSA +LT+ Cephase – ICR	with neutralizers lecithin, Tween® 80, β-lacatamases for inactivation of penicillins, all generations of cephalosporins and carbapenems	20 120	1.46076.0020 1.46076.0120
TSA with Penase LTHThio sedi. – ICR	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate, penase for inactivation of penicillins	20 120	1.46799.0020 1.46799.0120
ICR Settle Plates	150 mm , triple-bagged, gamma-irradiated, non-lockable		
TSA LTHThio sedi. – ICR 150 mm	3 plates, each triple-bagged – with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate	18	1.46820.0018



ICR/ICR+ Settle Plates

for isolators and critical cleanrooms



Description

Benefits

Related Products

Further Reading

Ordering Information

Product	Description	Content	Ord. No.
ICRplus Settle Plates	90 mm , triple-bagged, gamma-irradiated, lockable		
SDA +LTHTh – ICR+	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate	20	1.46702.0020
		120	1.46702.0120
TSA +LTHTh 90 mm ICR+	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate	20	1.46683.0020
		120	1.46683.0120
TSA + LTHThio sedi. – ICR+	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate	20	1.46787.0020
		120	1.46787.0120
TSA +LT 90 mm ICR+	with neutralizers lecithin, Tween® 80	20	1.46684.0020
		120	1.46684.0120
TSA 90 mm ICR+		20	1.46685.0020
		120	1.46685.0120
TSA +LT +Cephase – ICR+	with neutralizers lecithin, Tween® 80	120	1.46700.0120



ICR/ICR+ Settle Plates

for isolators and critical cleanrooms



Description

Benefits

Related Products

Further Reading

Ordering Information

Product	Description	Content	Ord. No.
ICRplus Settle Plates	90 mm , triple-bagged, gamma-irradiated, lockable		
TSA with Penase LTHThio sedi. – ICR+	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate	20	1.46800.0020
		120	1.46800.0120
Veg. Pept. Agar with LTHThio sedi. – ICR+	with neutralizers lecithin, Tween® 80, histidine, sodium thiosulfate	20	1.46804.0020
		120	1.46804.0120



LI Settle Plates

for microbial monitoring of air, surfaces and personnel



Description

Benefits

Related Products

Ordering Information



Microbial monitoring of air, surfaces and personnel requires culture media that ensure reproducible detection results. RT tryptic soy agar (TSA) and sabouraud dextrose agar (SDA) reliably detect bacteria, yeasts and molds in non-critical areas.

To monitor efficiently the status of environmental control, the culture media must be optimized for recovery of low levels of environmental contaminants.

LI settle plates contain a volume of 30 mL to reduce desiccation during sampling and incubation. Tryptic Soy Agar (TSA) – LI and Sabouraud Dextrose Agar (SDA) – LI are especially designed for active as well as passive air sampling. **LI settle plates** are optimized for long incubation periods.

[Link to our Website](#)



[Download Brochure](#)



[Request a demo](#)



LI Settle Plates



for microbial monitoring of air, surfaces and personnel

Description

Benefits

Related Products

Ordering Information

- Controlled high-quality raw materials
- **Data Matrix barcode** to ensure full traceability of each individual plate
- Filled under aseptic conditions
- Long shelf life to reduce incoming goods controls
- Storage at room temperature to keep at site of use
- SDA in colored plates to minimize risk of mix-ups
- High filling volume to allow prolonged exposure and incubation



LI Settle Plates

for microbial monitoring of air, surfaces and personnel



Description

Benefits

Related Products

Ordering Information

MAS air sampler family

The industry standard for viable air sampling.



MAS-100 NT®



MAS-100 VF®



MAS-100 Iso MH®



MAS-100 Iso NT®

Link to our Website



LI Settle Plates

for microbial monitoring of air, surfaces and personnel



Description

Benefits

Related Products

Ordering Information

Product	Description	Content	Ord. No.
LI Settle Plates	90 mm , single-bagged, non-lockable		
SDA – LI acc. EP		20	1.46028.0020
		120	1.46028.0120
SDA + Chloramphenicol – LI 30 mL		20	1.46003.0020
		120	1.46003.0120
SDA +LTHTTh – LI 30 mL	with neutralizers lecithin, Tween [®] 80, histidine, sodium thiosulfate	20	1.46052.0020
		120	1.46052.0120
TSA – LI 30 mL EP+USP		20	1.46004.0020
		120	1.46004.0120
TSA +LTHTTh – LI 30 mL	with neutralizers lecithin, Tween [®] 80, histidine, sodium thiosulfate	20	1.46002.0020
		120	1.46002.0120



ACTIVE VIABLE AIR MONITORING



MAS-100 NT®

The industry standard for viable air sampling

Cleanroom Grade
A/B

Cleanroom Grade
C/D



MAS-100 VF®

The easy to use air sampler

Cleanroom Grade
A/B

Cleanroom Grade
C/D



MAS-100 Iso MH®

Pure innovation

Isolator / RABS



MAS-100 Iso NT®

For isolators

Isolator / RABS

MAS-100 NT[®]

The industry standard for viable air sampling



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information



The **MAS-100 NT[®]** is the world's best-selling high performance air sampler and the premium version of the portable and battery-operated MAS-100[®] instruments. It is used for microbial monitoring of air and stands for highest safety and sampling efficiency.

- Integrated mass flow sensor with temperature and pressure compensation regulate the flow rate to 100 SLPM
- Several models (e.g., with HEPA exhaust filter and/or devices suitable for explosion hazard areas) are available.

Link to our Website



Download Brochure



Request a demo



MAS-100 NT[®]

The industry standard for viable air sampling



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information



Swiveling head for horizontal unidirectional air flows

- Independently validated according to ISO 14698 (complies with EN 17141)
- Built-in air flow compensation guarantees precise and reproducible flow rates
- Suitable for ISO 5/GMP Grade A cleanrooms
- Minimized turbulences in unidirectional air flows
- Compatible with 90 mm Petri dishes, 55 mm contact plates and other plate formats
- Automated calibration traceable to ISO 17025 / NIST
- Suitable for disinfection with vaporized and gaseous hydrogen peroxide and many industry-standard disinfectant

MAS-100 NT[®]



The industry standard for viable air sampling

Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

General

Operation principle	Sieve impaction (Andersen Sampler)
Flow rate	Default: 100 SLPM \pm 2.5%
Flow regulation	In-built mass flow sensor
Sample volumes	Selectable from 10 to 2000 liters Preset: 50, 100, 250, 500 & 1000 liters
Sequential Sampling Mode (Intervals)	40 to 2000 liters in 2-50 fractions up to 23:59 hours
Delay function	1 to 120 min
Nominal d_{50} value	1.1 μ m at 100 SLPM (confirmed by independent validation study acc. ISO 14698)
Biological efficiency	\sim 83% (independent validation study acc. ISO 14698)
Particle generation	Compliant with ISO 14644 1, Class 5



MAS-100 NT[®]



The industry standard for viable air sampling

Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Housing

Material	Anodized aluminum
Dimensions (H x W x L)	270 x 109 x 170 mm (with handle)
Weight	2.37 kg (with perforated lid)
Tripod screw thread	3/8" UNC (and 1/4" UNC insert), for use with standard camera tripod

Perforated lid

Material	Anodized aluminum
Dimensions (D x H)	Ø 109 x 32 mm
Weight	190 g
Lid geometry	Default: 300 x 0.6 mm
Sterilization	Autoclaving; dry heat

Others

Battery capacity	Approx. 5.5h or 33.000 liters (2 h or 12.000 liters for MAS-100 NT Ex [®] versions)
ATEX certification	II 3G Ex ec IIC T4 Gc (for MAS-100 NT Ex [®] version only)



MAS-100 NT[®]



The industry standard for viable air sampling

Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

ICR/ICRplus Contact Plates

Gamma-irradiated and triple-bagged contact plates are suitable for use in critical cleanrooms, RABS and isolators.



[Link to our Website](#)



MAS-100 NT[®]

The industry standard for viable air sampling



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Product	Description	Ord. No.
Perforated lid	300 holes x 0.6 mm, aluminum	1.09195.0001
Perforated lid	300 holes x 0.6 mm, aluminum, edge protection (MTO)*	1.19363.0001
Dust cover	Aluminum	1.09084.0001
Battery charger with Lemo plug	For all MAS-100 NT Ex [®] and MAS-100 NT [®] with serial numbers < 110499	1.09200.0001
Battery charger with DC plug	For MAS-100 NT [®] with serial number ≥ 110500	1.19133.0001
Tripod	For sampling in heights (88 - 250 cm)	1.09326.0001
Quick Change tripod adapter	For quick adaption of MAS-100 NT [®] to tripods without screwing	1.09223.0001
Mounting kit for filter	For upgrading MAS-100 NT [®] /Ex [®] with HEPA exhaust filter	1.17276.0001
HEPA H13 replacement filter	Ø 74 mm; Filter efficiency H13	1.17278.0001
Tube adapter	For sampling hardly to access sampling locations, can be used with tubing	1.09224.0001



MAS-100 NT®



The industry standard for viable air sampling

Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Product	Description	Ord. No.
Perforated lid contact plates	300 holes x 0.6 mm, aluminum, for use with 55 mm ICR contact plates	1.19149.0001
Contact plate adapter	For centering and fixing 55 mm contact plates	1.09214.0001
RapidMicroBio perforated lid with clamps	300 holes x 0.6 mm, aluminum with clamps to center GD plates (MTO)*	1.17316.0001
RapidMicroBio perforated lid	300 holes x 0.6 mm, aluminum for GD plates (MTO)*	1.17315.0001
Perforated lid	400 holes x 0.7 mm, aluminum	1.09088.0001
Perforated lid	400 holes x 0.7 mm, aluminum, edge protected (MTO)*	1.19364.0001
Perforated lid contact plates	400 holes x 0.7 mm, aluminum, for use with 55 mm ICR contact plates (non-lockable only)	1.09213.0001
RapidMicroBio perforated lid with clamps	400 holes x 0.7 mm, aluminum with clamps for centering GD plates (MTO)*	1.19166.0001
RapidMicroBio perforated lid	400 holes x 0.7 mm, aluminum for GD plates (MTO)*	1.17318.0001

*MTO products are not in stock – please mind lead times



MAS-100 NT®



The industry standard for viable air sampling

Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Product	Description	Ord. No.
MAS-100 NT®	Standard microbial air sampler	1.09191.0001
MAS-100 NT® with Filter	MAS-100 NT® (built-in H13 filter)	1.17274.0001
MAS-100 NT Ex®	For explosion hazard areas	1.09194.0001
MAS-100 NT Ex® with Filter	For explosion hazard areas (built-in H13 filter) (MTO)*	1.17275.0001
Scope of delivery	MAS-100 NT® microbial air sampler including: 1 x perforated lid (300 x 0.6 mm) 1 x dust cover 1 x Allan key No. 3 1 x USB cable 1.5 m 1 x battery charger with DC plug and regional plug adapters (US/JPN, EU, GB, AUS) 1 x USB stick with manuals and Feller correction table calibration certificates for MAS-100 NT®/NT Ex® and calibration device 1 x transportation case	

*MTO products are not in stock – please mind lead times

MAS-100 VF[®]

The easy to use air sampler



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information



The **MAS-100 VF[®]** is a portable and battery-operated microbial air sampler. It has the same collection efficiency as the premium version MAS-100 NT[®].

- Compact design
- EasyClean Plate holder
- Speed-based air volume flow control

Link to our Website



Download Brochure



Request a demo



MAS-100 VF[®]

The easy to use air sampler



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information



Single touch slide control

- Independently validated according to ISO 14698 (complies with EN 17141)
- Configurable and programmable via common web browsers
- Compatible with 90 mm Petri dishes, 55 mm contact plates and other plate formats
- Calibration certificate stored in the instrument
- Automated calibration traceable to ISO 17025 / NIST

Learn more
on our Website



MAS-100 VF[®]

The easy to use air sampler



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

General

Operation principle	Sieve Impaction
Flow rate	Default: 100 LPM \pm 4%
Flow regulation	Blower Speed Regulation
Sample volumes	Selectable from 10 to 2000 liters Preset: 50, 100, 250, 500 & 1000 liters
Delay function	1 to 120 min
Nominal d_{50} value	1.1 μ m at 100 LPM (confirmed by independent validation study acc. ISO 14698)
Biological efficiency	\sim 77% (independent validation study acc. ISO 14698)

Housing

Material	Anodized aluminum
Dimensions (H x W x L)	179 x 109 x 148 mm (with handle)
Weight	1.92 kg (with perforated lid)
Tripod screw thread	$\frac{3}{8}$ " UNC (and $\frac{1}{4}$ " UNC insert), for use with standard camera tripod



MAS-100 VF[®]

The easy to use air sampler



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Sampling Head (Perforated lid)

Material Anodized aluminum

Dimensions (D x H) Ø 109 x 32 mm

Weight 190 g

Lid geometry Default: 300 x 0.6 mm

Sterilization Autoclaving; dry heat

Others

Battery capacity Approx. 5 h or 30.000 liters



MAS-100 VF[®]

The easy to use air sampler



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

ICR/ICRplus Settle Plates

Our gamma-irradiated and triple-bagged settle plates for monitoring of isolators and critical cleanrooms.



Link to our Website



MAS-100 VF[®]

The easy to use air sampler



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Product	Description	Ord. No.
Perforated lid	Default version: 300 holes x 0.6 mm, aluminum	1.09195.0001
Perforated lid contact plates	300 holes x 0.6 mm, aluminum, for use with 55 mm ICR contact plates	1.19149.0001
Dust cover	Aluminum cover protects sieve from dirt	1.09084.0001
Battery charger	Including regional adapters and USB cable (USB-C to Type B socket)	1.17104.0001
Tripod	For sampling in heights (88 - 250 cm)	1.09326.0001
Quick change tripod adapter	For quick adaption of MAS-100 VF [®] to tripods without screwing	1.09223.0001
Tube adapter	For sampling hardly to access sampling locations, can be used with tubing	1.09224.0001



MAS-100 VF®

The easy to use air sampler



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Product	Description	Ord. No.
Perforated lid	300 holes x 0.6 mm, aluminum, edge protection (MTO)*	1.19363.0001
RapidMicroBio perforated lid	300 holes x 0.6 mm, aluminum for GD plates (MTO)*	1.17315.0001
Perforated lid	400 holes x 0.7 mm, aluminum	1.09088.0001
Perforated lid	400 holes x 0.7 mm, aluminum, edge protected (MTO) *	1.19364.0001
Perforated lid contact plates	400 holes x 0.7 mm, aluminum, for use with 55 mm ICR contact plates (non-lockable only)	1.09213.0001
RapidMicroBio perforated lid	400 holes x 0.7 mm, aluminum for GD plates (MTO)*	1.17318.0001
Petri dish clamps (3x)	With screws	1.17171.0001

***MTO products are not in stock – please mind lead times**



MAS-100 VF®

The easy to use air sampler



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Product	Description	Ord. No.
MAS-100 VF®	Microbial air sampler	1.17103.0001
Scope of delivery	MAS-100 VF® microbial air sampler including: 1 x perforated lid (300 x 0.6 mm) 1 x dust cover (plastic) 1 x Allan key No. 3 1 x USB cable 1.5 m 1 x battery charger with DC plug and regional plug adapters 1 x USB stick with manuals and Feller correction table calibration certificates for MAS-100 VF® and calibration device 1 x transportation case	

MAS-100 Iso MH[®]

Pure innovation



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

The **MAS-100 Iso MH[®]** is a high-performance, remote microbial air sampler designed for use with isolators and RABS systems. Its sampling principle, flow rate regulation using a mass flow sensor, accuracy and collection efficiencies are identical to MAS-100 NT[®]. All instruments are delivered with standard interfaces (USB, RS232, SPS and 9 digital I/O).

- Up to four fixed sample locations within the isolator
- Up to 10 m (50 m without decontamination) distance between sampling location inside and control unit
- Models with additional interfaces such as Ethernet, ProfiBus or ProfiNet are optional



Link to our Website



Download Brochure



Request a demo



MAS-100 Iso MH[®]

Pure innovation



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information



Unique double valve system

- Double valve system prevents back-flow and cross-contamination and allows active VHP decontamination
- High collection efficiency, independent validation according to ISO 14698 (equivalent to EN 17141)
- Built-in air flow compensation for precise, reproducible flow rate of 100 SLPM
- Various stainless steel sampling head configurations for 90 mm plates or other plate formats
- Multiple options for combining filters (in-line/exhaust)
- Automated calibration with MAS-100 Regulus traceable to ISO 17025 or NIST

MAS-100 Iso MH[®]

Pure innovation



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

General

Operation principle	Sieve Impaction (Andersen Sampler)
Flow rate	Default: 100 SLPM \pm 2.5%
Flow regulation	In-built mass flow sensor
Sample volumes	Selectable from 10 to 2000 liters Preset: 50, 100, 250, 500 & 1000 liters
SQS (sequential sampling mode, Intervals)	40 to 2000 liters in 2-50 fractions up to 23:59 hours
CNS Mode (consecutive sampling)	Order of sampling/decontamination automatically one head after the other
Delay function	1 to 120 min
Nominal d_{50} value	1.1 μ m at 100 SLPM (confirmed by independent validation study acc. ISO 14698)
Biological efficiency	\sim 78% (independent validation study acc. ISO 14698)



MAS-100 Iso MH[®]

Pure innovation



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Housing

Material	Side walls: anodized aluminum; body: electro polished stainless steel
Dimensions (H x W x L)	160 x 230 x 283 mm
Weight	7.5 kg

Sampling Head with perforated lid

Material	Stainless Steel
Dimensions (D x H)	Ø 109 x 90 mm
Weight	1.5 kg
Lid geometry	Default: 300 x 0.6 mm
Sterilization	Autoclaving; dry heat

Power

Input external power supply	110-240 V AC, 0.7-1.5 A / 47-63 Hz
Input instrument	24 V DC ± 5%, Max. Limits: 2.7 A and 65 W



MAS-100 Iso MH[®]

Pure innovation



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information



IsoBag[®] rapid transfer bag

A unique transportation and packaging bag as a reliable and easy way to have environmental monitoring plates right at hand in the isolator.

Link to our Website



MAS-100 Iso MH[®]

Pure innovation



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Product	Description	Ord. No.
Standard base for sampling Head	Stainless steel with Allen key, tri-clamp, and silicone gasket	1.09328.0001
Silicone gasket (5x)	Sealing between standard base for sampling head and tubing	1.17084.0001
Easy clean base for sampling head	Stainless steel with tri-clamp and silicone screen gasket	1.17091.0001
Silicone screen gasket (5x)	Sealing between easy clean base for sampling head and tubing	1.17099.0001
Perforated lid	Default version: 300 holes x 0.6 mm, edge-protected, stainless steel	1.09189.0001
Dust cover	Stainless steel cover protects sieve from dirt	1.09644.0001
Power supply	With power cord and regional plug adapters	1.17182.0001
MAS-100 Iso NT [®] IP54 Kit	protection against solid objects that are bigger than 1 mm and water splashing from any direction	1.17159.0001
Remote cable MAS-100 Iso NT [®]	Can be used for remote start & stop of sampling and can display the instrument status via LED	1.17181.0001

Please mind lead times for accessories



MAS-100 Iso MH[®]

Pure innovation



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Product	Description	Ord. No.
Perforated lid	400 holes x 0.7 mm, edge-protected, stainless steel	1.09222.0001
RapidMicroBio easy clean base	Stainless steel, for Growth Direct EM Cassettes	1.17319.0001
RapidMicroBio perforated lid	300 holes x 0.6 mm, stainless steel for Growth Direct Cassettes	1.17314.0001
Aervent KTGR Opticap XL5 capsule	Hydrophobic 0.2 µm, 5''	KTGRA05TT1
Aervent KTGR Opticap XL10 capsule	Hydrophobic 0.2 µm, 10''	KTGRA10TT1
Pressure test kit	For testing tightness of installation	1.17085.0001
Calibration cable, 10 m	RS232 connection for connection of calibration device MAS-100 [®] Regulus	1.17098.0001
ProfiNet module	Converts ProfiBus to ProfiNet interface	1.17290.0001
Elbow joint	Stainless steel with 2 tri-clamps and silicone gaskets	1.17083.0001
Tri-clamp (3x)	To facilitate mounting of sampling heads	1.09440.0001

Please mind lead times for accessories



MAS-100 Iso MH[®]

Pure innovation



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Product	Interface	No. of heads	Ord. No.	Interface	No. of heads	Ord. No.
MAS-100 Iso MH[®] microbial air sampler	Standard*	1	1.17174.0001	ProfiBus*	1	1.17177.0001
		2	1.17118.0001		2	1.17144.0001
		3	1.17146.0001		3	1.17147.0001
		4	1.17149.0001		4	1.17155.0001
	Ethernet*	1	1.17178.0001	ProfiNet*	1	1.17292.0001
		2	1.17145.0001		2	1.17309.0001
		3	1.17148.0001		3	1.17310.0001
		4	1.17157.0001		4	1.17312.0001

Scope of delivery

MAS-100 Iso NT[®] microbial air sampler including:

MAS-100 Iso NT[®] control unit: The control units come with pre-assembled valves to connect optional tubing, sampling base and perforated lids.

*MTO products are not in stock – please mind lead times

MAS-100 Iso NT®

For isolators



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

The **MAS-100 Iso NT®** is a high-performance, remote microbial air sampler designed for use with isolators and RABS systems. Its sampling principle, flow rate regulation using a mass flow sensor, accuracy and collection efficiencies are identical to MAS-100 NT®. All instruments are delivered with standard interfaces (USB, RS232, SPS and 9 digital I/O).

- Alphanumeric LCD display
- One fixed sample locations within the isolator
- Up to 10 m (50 m without using decontamination) distance between sampling location inside and control unit
- Models with additional interfaces such as Ethernet, ProfiBus or ProfiNet are optional



Link to our Website



Download Brochure



Request a demo



MAS-100 Iso NT[®]

For isolators



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information



Unique double valve system

- Double valve system prevents back-flow and cross-contamination and allows active VHP decontamination
- High collection efficiency, independent validation according to ISO 14698 (equivalent to EN 17141)
- Built-in air flow compensation for precise, reproducible flow rate of 100 SLPM
- Various stainless steel sampling head configurations for 90 mm plates or other plate formats
- Multiple options for combining filters (in-line/exhaust)
- Automated calibration with MAS-100 Regulus traceable to ISO 17025 or NIST

MAS-100 Iso NT[®]

For isolators



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

General

Operation principle	Sieve Impaction (Andersen Sampler)
Flow rate	Default: 100 SLPM \pm 2.5%
Flow regulation	In-built mass flow sensor
Sample volumes	Selectable from 10 to 2000 liters Preset: 50, 100, 250, 500 & 1000 liters
SQS (sequential sampling mode, Intervals)	40 to 2000 liters in 2-50 fractions up to 23:59 hours
Delay function	1 to 120 min
Nominal d_{50} value	1.1 μ m at 100 SLPM (confirmed by independent validation study acc. ISO 14698)
Biological efficiency	\sim 78% (independent validation study acc. ISO 14698)



MAS-100 Iso NT[®]

For isolators



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Housing

Material	Side walls: anodized aluminum; body: electro polished stainless steel
Dimensions (H x W x L)	160 x 230 x 283 mm
Weight	7.5 kg
Display	Alphanumeric LCD

Sampling Head with perforated lid

Material	Stainless Steel
Dimensions (D x H)	Ø 109 x 90 mm
Weight	1.5 kg
Lid geometry	Default: 300 x 0.6 mm
Sterilization	Autoclaving; dry heat

Power

Input external power supply	110-240 V AC, 0.7-1.5 A / 47-63 Hz
Input instrument	24 V DC ± 5%, Max. Limits: 2.7 A and 65 W



MAS-100 Iso NT[®]

For isolators



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information



IsoBag[®] rapid transfer bag

A unique transportation and packaging bag as a reliable and easy way to have environmental monitoring plates right at hand in the isolator.

Link to our Website



MAS-100 Iso NT[®]

For isolators



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Product	Description	Ord. No.
Standard base for sampling Head	Stainless steel with Allen key, tri-clamp, and silicone gasket	1.09328.0001
Silicone gasket (5x)	Sealing between standard base for sampling head and tubing	1.17084.0001
Easy clean base for sampling head	Stainless steel with tri-clamp and silicone screen gasket	1.17091.0001
Silicone screen gasket (5x)	Sealing between easy clean base for sampling head and tubing	1.17099.0001
Perforated lid	Default version: 300 holes x 0.6 mm, edge-protected, stainless steel	1.09189.0001
Dust cover	Stainless steel cover protects sieve from dirt	1.09644.0001
Power supply	With power cord and regional plug adapters	1.17182.0001
MAS-100 Iso NT [®] IP54 Kit	protection against solid objects that are bigger than 1 mm and water splashing from any direction	1.17159.0001
Remote cable MAS-100 Iso NT [®]	Can be used for remote start & stop of sampling and can display the instrument status via LED	1.17181.0001

Please mind lead times for accessories



MAS-100 Iso NT[®]

For isolators



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Product	Description	Ord. No.
Perforated lid	400 holes x 0.7 mm, edge-protected, stainless steel	1.09222.0001
RapidMicroBio easy clean base	Stainless steel, for Growth Direct EM Cassettes	1.17319.0001
RapidMicroBio perforated lid	300 holes x 0.6 mm, stainless steel for Growth Direct Cassettes	1.17314.0001
Aervent KTGR Opticap XL5 capsule	Hydrophobic 0.2 µm, 5''	KTGRA05TT1
Aervent KTGR Opticap XL10 capsule	Hydrophobic 0.2 µm, 10''	KTGRA10TT1
Pressure test kit	For testing tightness of installation	1.17085.0001
Calibration cable, 10 m	RS232 connection for connection of calibration device MAS-100 [®] Regulus	1.17098.0001
ProfiNet module	Converts ProfiBus to ProfiNet interface	1.17290.0001
Elbow joint	Stainless steel with 2 tri-clamps and silicone gaskets	1.17083.0001
Tri-clamp (3x)	To facilitate mounting of sampling heads	1.09440.0001

Please mind lead times for accessories



MAS-100 Iso NT[®]

For isolators



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Product	Description	Interface	No. of heads	Ord. No.
MAS-100 Iso NT[®]	Microbial air sampler	Standard	1	1.09168.0001
		Ethernet	1	1.09174.0001
		ProfiBus*	1	1.09173.0001
		ProfiNet*	1	1.17291.0001

Scope of delivery

MAS-100 Iso NT[®] microbial air sampler including:

MAS-100 Iso NT[®] control unit: The control units come with pre-assembled valves to connect optional tubing, sampling base and perforated lids.

*MTO products are not in stock – please mind lead times



COMPRESSED GAS MONITORING



MAS-100 ATMOS®
The NEW microbial
compressed gas sampler

**Cleanroom Grade
A/B**

**Cleanroom Grade
C/D**



MAS-100 ATMOS®

The NEW microbial compressed gas sampler



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

The **MAS-100 Atmos**® microbial compressed gas sampler is a high precision, portable and battery-powered device to monitor pressurized gases such as air, nitrogen, carbon dioxide, and argon. The instrument is easy to use and can be operated in a pressure range of 1.2 to 7 bar.a (17.4 to 101.5 psi).

At flow rate of 100 liter the **MAS-100 Atmos**® has the identical impaction speed and D50 value as the MAS-100 NT®.

The microorganisms are impacted under pressure to avoid any damage by dynamic decompression.



Link to our Website



Request a demo



MAS-100 ATMOS®

The NEW microbial compressed gas sampler



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information



**Impacts microorganisms
under pressure**

Ease of use that minimizes handling errors

- Automated flow rate regulation for pre-defined gas types
- Automated decompression
- Automated calibration traceable to ISO 17025 / NIST
- Impaction under pressure minimizes risk of false negatives

From non-digital to fully computerized acc. to 21 CFR Part 11

- Freely configurable out-of-the-box use
- Individual login via hardware key token and PIN
- Tamper-proof audit trail download
- Unique Result barcode for implementation in EM software systems

MAS-100 ATMOS[®]

The NEW microbial compressed gas sampler



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

General

Operation principle	Sieve Impaction
Gas types	Default: compressed air, Nitrogen (N ₂), carbon dioxide (CO ₂), Argon (Ar), customized gas type (non-oxidizing or flammable gases) if gas factor known
Pressure range	1.2 to 7 bar.a ± 2% (17.4 to 101 psi)
Flow rate	Default: 100 LPM ± 5% Optional: 50 LPM ± 5%
Flow regulation	Mass-flow sensor
Sample volumes	Selectable from 50 to 3000 normal liters (NL) Preset: 50, 100, 250, 500 & 1000 NL
Nominal d ₅₀ value	1.1 µm at 100 LPM

Housing

Material	Injection-molded high-performance polystyrene with plastic primer and texture varnish
Dimensions (H x W x L)	180 x 210 x 500 mm (with handle)
Weight	7.3 kg (with sampling head)



MAS-100 ATMOS®



The NEW microbial compressed gas sampler

Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Perforated lid

Material	Anodized aluminum (blue)
Dimensions (D x H)	Ø 109 x 90 mm
Weight	190 g
Lid geometry	Default: 300 x 0.6 mm
Sterilization	Autoclaving; dry heat

Others

Battery capacity	Approx. 2 to 5h depending on chosen settings and pressure range
Operating temperature	+5 to +45°C
Sound level	< 80 dB at highest pressure point (7 bar.a)



MAS-100 ATMOS®

The NEW microbial compressed gas sampler



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

ICR/ICRplus Settle Plates

Our gamma-irradiated and triple-bagged settle plates for monitoring of isolators and critical cleanrooms.



Link to our Website



MAS-100 ATMOS®

The NEW microbial compressed gas sampler



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Product	Description	Ord. No.
Sampling head Atmos with perforated lid 300 x 0.6 mm	Default version: provides D ₅₀ of 1.1 µm at 100 liters/min	1.17357.0001
Sampling head Atmos with perforated lid 300 x 0.47 mm	Provides D ₅₀ of 1.1 µm at 50 liters/min	1.17363.0001
Pressure tube Atmos 2 m with TC Connector	PTFE Tube with 3/8'' female connector to gas supply with connector to instrument	1.17354.0001
Gas exhaust tube Atmos 5 m with TC Connector	PVC tube Ø 50 mm, stainless steel connector and tri-clamp	1.17349.0001
Adapter gas exhaust tube for MAS-100 Atmos®	Stainless steel adapter for self-assembly of gas exhaust tubing	1.17348.0001
Hardware key set (5x) for MAS-100 Atmos®	Re-writable for individual user login	1.17333.0001
Wheeled transport case for MAS-100 Atmos®	For safe transportation	1.17334.0001

Please mind lead times for accessories



MAS-100 ATMOS®

The NEW microbial compressed gas sampler



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Product	Description	Ord. No.
Dust cover for MAS-100 Atmos®	Polypropylene, white, autoclavable	1.17340.0001
Plate holder for MAS-100 Atmos® stainless steel	Autoclavable	1.17335.0001
Filter cover for MAS-100 Atmos® with 3 screws	Aluminum, autoclavable	1.17336.0001
Tri-clamp 34 mm for MAS-100 Atmos®	Used to connect pressure tube with sampling head	1.17367.0001
Silicone gasket (5x) 34/10 mm for MAS-100 Atmos®	Sealing of connection between pressure gas tube and sampling head	1.17364.0001
Silicone gasket 50.5/40 mm for MAS-100 Atmos®	Sealing of connection between gas exhaust adapter and tubing	1.17366.0001
Battery charger MAS-100 Atmos®	Including USB-C cable and regional plug adapters	1.17329.0001

Please mind lead times for accessories



MAS-100 ATMOS®

The NEW microbial compressed gas sampler



Description

Benefits

Technical specifications

Related Products

Accessories

Ordering Information

Product	Description	Ord. No.
MAS-100 Atmos®	Microbial compressed gas sampler	1.17328.0001
Scope of delivery	MAS-100 Atmos® microbial compressed gas sampler including: 1 x sampling head with perforated lid (300 x 0.6 mm) 1 x pressure connection tube 2 m PTFE with connection pieces (3/8" female) to gas supply and instrument (34 mm tri-clamp and silicone gasket) 1 x dust lid 1 x tri-clamp Ø 34 mm 1 x silicone gasket 34/10 mm 1 x battery charger incl. USB C cable and regional plug adapters 1 x system administrator hardware key (black) 5 x configurable hardware keys documentation incl. calibration certificates	



MEDIA FILL / ASEPTIC PROCESS SIMULATION



**Dehydrated
Culture Media (DCM)**

**Cleanroom Grade
A/B**

Isolator / RABS



**Ready-to-use
media fill**

**Cleanroom Grade
A/B**

Isolator / RABS



Dehydrated Culture Media (DCM)

for safe media preparation



Description

Benefits

Related Products

Ordering Information

The available **granulated culture media** for media fill tests are irradiated and triple-packaged. Every batch is carefully tested to ensure excellent growth properties and filterability, and of course absence of microbial contamination.

The granulation format of the media reduces dust spread during preparation, thus lowering health risks and preventing contamination of the work environment. Less dust on equipment also eases handling and cleaning.

Tryptic soy broth and its non-animal origin alternative are available in 500 g and 5 kg packs.



[Link to our Website](#)



[Download Brochure](#)



Dehydrated Culture Media (DCM)

for safe media preparation



Description

Benefits

Related Products

Ordering Information

- **Safe media preparation** thanks to irradiated and triple-bagged products
- **Low-dust granules** to maximize safety and performance
- **Excellent cold filterability** to avoid clogging and costs
- **Excellent solubility** to significantly reduce preparation time and effort
- BSE-certified and non-animal origin products for **ultimate process safety**



Dehydrated Culture Media (DCM)

for safe media preparation



Description

Benefits

Related
Products

Ordering
Information

RTU bags

We also have Tryptic Soy Broth and Vegetable peptone broth available as ready to use media for media fill tests



Dehydrated Culture Media (DCM)



for safe media preparation

Description

Benefits

Related Products

Ordering Information

Product	Content	Ord. No.
Tryptic Soy Broth non-animal origin	500 g	1.00550.0500
	5 kg	1.00550.5000
Tryptic Soy Broth	500 g	1.00800.0500
	5 kg	1.00800.5000

Ready-to-use media fill

Pre-filtered and ready to go



Description

Benefits

Related Products

Ordering Information

10 L media fill bags containing ready-to-use culture media

The 10 L bags are filled with either tryptic soy broth or vegetable peptone broth. The media is sterile, triple-bagged and irradiated for safe use in process lines. The self-collapsing bags include an 80 cm long tube with a connector to enable direct connection to the filling line. The bags are available with an MPC or a Lynx connector.



Link to our Website



Download Brochure



Ready-to-use media fill

Pre-filtered and ready to go



Description

Benefits

Related Products

Ordering Information

- **Ready and easy to use** with simple connection to your process line
- Packed in transport box for easy handling
- Sterile, triple-bagged and irradiated for **ultimate process safety**
- BSE-certified and non-animal origin products for **total confidence**



Ready-to-use media fill

Pre-filtered and ready to go



Description

Benefits

Related Products

Ordering Information

DCM

Our media for media fill testing is as well available in granulated and irradiated form



Ready-to-use media fill



Pre-filtered and ready to go

Description

Benefits

Related Products

Ordering Information

Product	Content	Ord. No.
Tryptic Soy Broth	10 L bag	1.46316.0001
Vegetable Peptone Broth	10 L bag	1.46332.0001
Vegetable Peptone Broth with Lynx connector	10 L bag	1.46749.0001

A/B

CLEANROOM GRADE



**ICR/ICR+
Contact Plates**



**ICR
swabs**



**HYCON®
Contact Slides**



**ICR/ICR+
Settle Plates**



**IsoBag® Rapid
Transfer Bag**



**MAS-100
NT®**



**MAS-100
Atmos®**



**MAS-100
VF®**



**Dehydrated
Culture Media**



**Ready-to-use
media fill**



CLEANROOM GRADE



**RT/RT+
Contact Plates**



**LI
Settle Plates**



**HYCON®
Contact Slides**



**MAS-100
NT®**



**MAS-100
Atmos®**



**MAS-100
VF®**



ISOLATOR / RABS



**ICR/ICR+
Contact Plates**



**ICR/ICR+
Settle Plates**



**IsoBag® Rapid
Transfer Bag**



ICR swabs



**MAS-100
Iso NT®**



**MAS-100
Iso MH®**



**Dehydrated
Culture Media**



**Ready-to-use
media fill**



CONTACT



**Your first choice for
Environmental Monitoring Solutions
for Pharmaceutical Industry:**

**Learn more on
our Website**



**For further information
please contact:**

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany

[SigmaAldrich.com](https://www.sigmaaldrich.com)

© 2022 Merck KGaA, Darmstadt, Germany and its affiliates. All Rights Reserved.
Merck, the vibrant M, Millipore, IsoBag, Opticap and HYCON are trademarks of Merck KGaA, Darmstadt,
Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information
on trademarks is available via publicly accessible resources.

Lit. No. MK_iPDF9514EN Ver. 1.0
07/2022